

LIGHTWEIGHT PORTABLE SCOOTER



Seat Color  
Black Vinyl

Body Colors  
Blue  
Red



3 wheel GB117D  
4 wheel GB147D

The next generation of the popular Buzzaround scooter series is now available! The Buzzaround XL offers more standard features to make your riding experience better than ever. Features a wraparound Delta tiller, a new LED battery gauge and a new LED angle-adjustable headlight, the new Buzzaround XL is also longer to offer more foot room, with a great turning radius. The Buzzaround XL disassembles in seconds with just four easy steps making transporting it from one destination to another a breeze!





## Wireless disassembly in 4 easy steps!

Item	GB117D	GB147D
Description	XL	XL
Medicare Code	K0800	K0800
Weight Capacity	300 lbs.	300 lbs.
Drive Wheels	Rear	Rear
Maximum Speed	4 mph	4 mph
Operating Range up to <sup>1</sup>	8 miles	8 miles
Under Deck Clearance	3"	3"
Turning Radius	37"	47"
Type Batteries	2-12AH	2-12AH
Colors	Red, Blue	Red, Blue
Freewheel Mode	Yes	Yes
Electronic Speed Control	Yes	Yes
Electro-Mechanical Brakes	Yes	Yes
Charger	Off-board	Off-board
Controller	50 Amp Dynamic R-Series	50 Amp Dynamic R-Series

Dimensions		
Length	42.5"	41.5"
Width	21.5"	21.5"
Height (ground to top of back of seat)	33"-35"	33"-35"
Ground to Top of Seat	21"-23"	21"-23"
Ground to Top of Deck	4.75"	4.75"
Deck to Top of Seat	16"-18"	16"-18"
Number of Seat Height Positions	3	3
Size of Increments	1"	1"
Front Axle to Rear Axle	33.5"	32"

Weight of Unit: (assembled) <sup>2</sup>		
Front Half	101 lbs.	109 lbs.
Rear Half	35 lbs.	43 lbs.
Seat w/Arms	25 lbs.	25 lbs.
Batteries	22 lbs.	22 lbs.
	19 lbs. pack	19 lbs. pack

Tires:		
Front	Solid	Solid
Rear	8"	8"
Rear Anti-Tip Wheels	8"	8"
	2"	2"

Standard Seat		
Seat Back Height - no headrest	Stadium Style	Stadium Style
Width x Depth	16"	16"
Color	18" x 16"	18" x 16"
	Black Vinyl	Black Vinyl
Optional Seat:		
Back Height - no headrest	Stadium Style	Stadium Style
Width x Depth	16"	16"
Color	20" x 18"	20" x 18"
	Black Vinyl	Black Vinyl

Warranty:		
Frame	Lifetime Warranty	Lifetime Warranty
Drive Train	2 Year Warranty	2 Year Warranty
Electronics <sup>2</sup>	2 Year Warranty	2 Year Warranty
Batteries	1 Year Warranty	1 Year Warranty

<sup>1</sup> Battery range tested at 200 lbs., and will vary due to rider weight, drive surface and terrain.

<sup>2</sup> Electronics warranty excludes batteries.



**Step 1:** Pop off the seat!



**Step 2:** Pull out the battery pack!



**Step 3:** Fold down tiller and lock in place!



**Step 4:** Pull forward on the drivetrain release lever to separate the front and rear frame sections!



Disassembles into 5 pieces!

4 wheel model shown



A battery charging port is conveniently located on the tiller under the control panel.

## Features

- 300 lb. weight capacity
- Tiller mounted basket
- Angle-adjustable LED headlight
- More Leg room and foot space