

## Ventura 4 DLX

HCPCS:




FIGURE A




FIGURE B



FIGURE C




FIGURE D



FIGURE E

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- 33AH batteries for extended range
- Ergonomic throttle control
- Height-adjustable swivel seat adjusted forward and backward features fold-down backrest and adjustable headrest
- Front and rear suspension with Anti-tip wheels
- Programmable controller
- Available with 18" or 20" width Captain's Seat
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust Delta tiller
- Headlight, rearview mirrors and turn signals/hazard lights
- Easy free wheel operation

SPECIFICATIONS	Description
Base Weight	99 lbs.
Batteries (pair)	12V x 33AH
Battery Weight (pair)	46 lbs.
Brakes	Electromagnetic
Charger	4A Offboard
Climbing Angle	7.5°
Dimensions	48" (L) x 25" (W)
Floor to Seat Height	22" - 24"
Front Wheels	10" x 3" Flat-free
Ground Clearance	2.5"
HCPCS	N/A
Heaviest Piece Weight	70 lbs.
Maximum Range	16 miles
Motor	24V x 350W x 5100 rpm
Rear Wheels	10" x 3" Flat-free
Seat Dimensions	18" (W) x 18" (D)   20" (W) x 18" (D)
Seat Type	Captain's
Seat Weight	40 lbs.
Top Speed	6 mph

SPECIFICATIONS	Description
Total Weight (excl. hand control)	185 lbs.
Turning Radius	54"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	400 lbs
Weight Capacity (13046   13046-HD)	400 lbs.

ACCESSORIES	Description
SF8030	Crutch Holder Set, SType, 1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8050	Sun Shade for Scooter,1/ea, EA
SF8020	Rear Basket Set, LType, 1/ea, EA
SF8030	Crutch Holder Set, SType, 1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8050	Sun Shade for Scooter,1/ea, EA
SF8020	Rear Basket Set, LType, 1/ea, EA

PARTS	Description
VENT-12	Rear Shroud, Folding Seat Ventura,1/ea, EA
VENT-13	StemBearingAssemblyVentura3&4,1/ea, EA
VENT-16	Anti-Tipper Wheel, Ventura 3 & 4,1/ea, EA
VENT-44	Head Light Assembly, Ventura 4, 1/ea, EA
VENT-17	FrontFrameVentura4,1/ea, EA
VENT-14	FrontBumperAssemblyVentura4,1/ea, EA
VENT-15	ShockAbsorberVentura3&4, EA
VENT-10	FrontShroudVentura4,1/ea, EA
VENT-11	MiddleShroudVentura3&4,1/ea, EA
VENT-08	ControllerVentura3&4,1/ea, EA
VENT-09	BatteryStrapsVentura,CPN,1/set, PR
VENT-06	FrontTillerCoverVentura3&4,1/ea, EA
VENT-07	RearTillerCoverVentura3&4,1/ea, EA
VENT-04	FacepltW/Switches18FshSeatVentura3&4/1ea, EA
VENT-05	FacepltW/Switches18/20CapSeatVentura3&4, EA
VENT-02	TillerSwivelAssemblyVentura3&4,1/ea, EA
VENT-03	HandlebarAssemblyVentura3&4,1/ea, EA
VENT-01	Forward/ReverseAssembly,Ventura3&4,1/ea, EA
VENT-37	WireHarnessVentura3&4,18"FishOnSeat,1/, EA
VENT-42	SeatAssemblyVentura3&4,18"CaptSeat,1/e, EA
VENT-12	Rear Shroud, Folding Seat Ventura,1/ea, EA
VENT-13	StemBearingAssemblyVentura3&4,1/ea, EA
VENT-16	Anti-Tipper Wheel, Ventura 3 & 4,1/ea, EA
VENT-44	Head Light Assembly, Ventura 4, 1/ea, EA
VENT-17	FrontFrameVentura4,1/ea, EA
VENT-14	FrontBumperAssemblyVentura4,1/ea, EA
VENT-15	ShockAbsorberVentura3&4, EA
VENT-10	FrontShroudVentura4,1/ea, EA
VENT-11	MiddleShroudVentura3&4,1/ea, EA
VENT-08	ControllerVentura3&4,1/ea, EA
VENT-09	BatteryStrapsVentura,CPN,1/set, PR
VENT-06	FrontTillerCoverVentura3&4,1/ea, EA

PARTS	Description
VENT-07	RearTillerCoverVentura3&4,1/ea, EA
VENT-04	FacepltW/Switches18FshSeatVentura3&4/1 ea, EA
VENT-05	FacepltW/Switches18/20CapSeatVentura3&4, EA
VENT-02	TillerSwivelAssemblyVentura3&4,1/ea, EA
VENT-03	HandlebarAssemblyVentura3&4,1/ea, EA
VENT-01	Forward/ReverseAssembly,Ventura3&4,1/ea, EA
VENT-37	WireHarnessVentura3&4,18"FishOnSeat,1/, EA
VENT-42	SeatAssemblyVentura3&4,18"CaptSeat,1/e, EA