

# ENHANCE

your Active Lifestyle  
with added Comfort  
and Performance



4-Wheel

**VICTORY® 10**



3-Wheel shown with  
optional 22" wide high-back seat

The **Victory® 10** series offers added style, comfort and performance for your active lifestyle. Standard features include LED pathway light, LED battery meter, feather touch disassembly, wraparound delta tiller, 10" solid tires and a front basket. A 400 lb. weight capacity, a per charge range of up to 16 miles and a maximum speed up to 5.3 mph compliment the versatility of the Victory 10. Available in Viper Blue and Candy Apple Red.



**PRIDE® MOBILITY SCOOTERS**  
LIVE YOUR BEST®

[pridemobility.com](http://pridemobility.com)





## Specifications 3W:SC610 4W:SC710

Weight Capacity	400 lbs.
Maximum Speed <sup>1</sup>	Up to 5.3 mph
Ground Clearance <sup>3</sup>	3W: 2.2" at motor, 1.26" at anti-tip 3.4" at side deck 4W: 2" at motor, 1.5" at anti-tip 3.1" at side deck
Turning Radius <sup>3</sup>	3W: 42" 4W: 54"
Overall Length <sup>3</sup>	3W: 45.25" 4W: 47"
Overall Width <sup>3</sup>	22.25"
Seat-to-Ground Height	20.5"-22"
Seat-to-Deck Height	17.5"-19"
Front Tires	3.5" x 10" solid
Rear Tires	3.5" x 10" solid
Suspension	N/A
Range Per Charge <sup>1,2</sup>	Up to 16 miles
Total Weight Without Batteries <sup>3,7</sup>	3W: 133 lbs. 4W: 137 lbs.
Heaviest Piece When Disassembled <sup>3</sup>	3W: 54 lbs. (front section) 4W: 61 lbs. (front section)
Standard Seating	Type: Foldable molded plastic Weight: 28 lbs. Material: Black Vinyl Dimensions: 18" width, 17" depth
Drive System	Rear-wheel drive sealed transaxle with 24-volt DC motor
Dual Braking System	Electronic, regenerative and electromechanical
Battery Requirements <sup>3,4,5</sup>	Size: (2) U-1 standard, Weight U-1: 24 lbs. each Size: (2) 40 AH (optional), Weight 40 AH: 32 lbs. each
Battery Charger	3.5-amp, off-board (U-1), 5-amp, off-board (40 AH)
Colors	 

Candy Apple Red Viper Blue

## HCPC Code: K0801 (Group1, Heavy Duty)

- <sup>1</sup> Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tire condition. This specification can be subject to a variance of (+ or -) 10%.
- <sup>2</sup> Tested in accordance with ANSI/RESNA, MC Vol2, Section 4 & ISO 7176-4 standards. Results derived from theoretical calculation based on battery specifications and drive system performance. Test conducted at maximum weight capacity.
- <sup>3</sup> Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- <sup>4</sup> AGM or gel cell type required.
- <sup>5</sup> Battery weight may vary based on battery manufacturer.
- <sup>6</sup> Includes standard seat weight.

## Features

- Exclusive one-hand feather touch disassembly provides the ultimate quick transport system
- Delta tiller with wraparound handles
- Easy-to-read LED battery meter
- Exclusive stylish, lightweight, non-scuffing, black, low-profile tires
- Long-lasting, high-intensity LED headlight offers optimal pathway illumination
- Lightweight pinchless hinge seat eight position rotation, sliders and viscoelastic foam inserts for increased comfort
- Easy-grip tiller adjustment knob
- Easy access tiller-mounted charger port

## Accessories

- Cane/crutch holder (*single, double*)
- Cell phone holder
- Cup holder
- Forearm crutch holder
- Oxygen tank holder
- Rear basket
- Saddlebag
- Safety flag
- Walker holder
- Weather cover
- XLR USB charger

## Options

- 40 AH batteries

## Optional Seating

- 18"W x 18"D black vinyl reclining high-back seat, with sliders, height adjustable headrest, flip-up, height and width adjustable armrest
- 22"W x 18"D black vinyl reclining high-back seat height adjustable headrest, flip-up, height and width adjustable armrest
- Optional seat sizes available

## Military Patches



## HCPC Code: K0801 (Group 1, Heavy Duty)



Pride offers a variety of lifts to **transport your mobility product with ease**. To find a compatible lift for your vehicle and mobility product, use our new Lift Finder application at [pridemobility.com/lift-finder/](http://pridemobility.com/lift-finder/)

NOTE: All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for scooter only. Overall measurements will vary based on seating and accessory selections.

(US) 800-800-8586

(CAN) 888-570-1113

[pridemobility.com](http://pridemobility.com)



**Pride**  
Mobility Products Corp.