

merits

EZ-GO

& EZ-GO DELUXE



EZ-GO Deluxe



EZ-GO

EZ-GO DELUXE

- Ergonomic high-back seat for additional back support: 18"W x 16"D x 18"H
- Higher weight capacity of 300 lbs.
- Quick connect feature for easy assembly and transportation
- Compact and lightweight design: 113 lbs. overall weight
- 22 Ah batteries for extended range
- Color options: pink, white and turquoise

EZ-GO

- Ergonomic high-back seat for additional back support: 18"W x 16"D x 18"H
- Quick connect feature for easy assembly and transportation
- Compact and lightweight design: 100 lbs. overall weight
- Turning radius of 23" for easy maneuverability in tight spaces
- Color options: red and blue

EZ-GO

HCPCS	K0813/K0814
Weight Capacity	250 lbs
Max. Speed	4 mph
Range	up to 9 miles
Turning Radius	23"
Ground Clearance	1.8"
Motor	DC 24V, 90W
Controller	PG VSI 50A
Batteries	12V/14Ah
Charger	1.5A off-board
Gradient	6°
Drive Wheel	8" x 3" solid
Casters	6" x 2"
Total Weight	100 lbs
Weight of Heaviest Piece	42 lbs

EZ-GO DELUXE

HCPCS	-
Weight Capacity	300 lbs
Max. Speed	4 mph
Range	up to 12 miles
Turning Radius	23"
Ground Clearance	2.25"
Motor	DC 24V, 90W
Controller	PG VSI 50A
Batteries	12V/22Ah
Charger	1.5A off-board
Gradient	6°
Drive Wheel	8" x 3" solid
Casters	7" x 1.75"
Total Weight	113 lbs
Weight of Heaviest Piece	54 lbs

OPTIONS

- Oxygen tank holder
- USB charger
- Cane and crutch holder
- Walker holder
- Foldable rear basket
- Cup holder



Merits Health Products Inc.

4245 Evans Avenue, Fort Myers FL, 33901

Phone: 1-800-963-7487 Fax: 239-772-32

merits