



TRAVELMATE™

THE WORLD IS YOURS TO EXPLORE



TRAVELMATE™



The TravelMate™ is designed to take with you wherever you go. It easily folds without disassembly, and fits in most modern vehicle trunks and cargo areas. As the most compact model in the Amigo lineup, TravelMate is perfect for people who want to keep moving while out and about.

- Travel compatibility:
- Travel surface: Indoor / outdoor hard flat surfaces
- Removable battery
- Compact, lightweight design



Twist-grip throttle
Features a twist-grip throttle that is intuitive and comfortable.

15"



Designed to go
The seat and handle fold down, allowing it to fit in most trunks and storage compartments.

— Specs —



250 lbs (113 kg)
Weight capacity

56 - 63 lbs (25.4 - 28.5 kg)
Unit weight (56 lbs. when battery removed)



24V 10Ah
Lithium-Ion battery



8 Miles (12.8 km)
On a full charge



0 - 4.5 MPH (0 - 7.24 kph)
Variable



Offboard charger
24V, cord kept separate from Amigo



38" x 22.25" (96.5 cm x 56.5 cm)
Unfolded dimensions



43" x 15" (109.22 cm x 38.1 cm)
Folded dimensions



3.75" (9.5 cm)
Ground clearance

— Color options —



Base blue

Base red

